

Clinical Insights from Navitor

Santiago Garcia, MD



TCT[®]

TRANSCATHETER
CARDIOVASCULAR
THERAPEUTICS[®]

Disclosure of Relevant Financial Relationships

Within the prior 24 months, I have had a financial relationship with a company producing, marketing, selling, re-selling, or distributing healthcare products used by or on patients:

Speaker Name: Santiago Garcia, MD

Nature of Financial Relationship

Grant/Research Support

Consultant Fees/Honoraria

Ineligible Company

Institutional Grants Edwards Lifesciences,
Abbott Structural Heart, Medtronic, BSCI

Edwards Lifesciences, Abbott Structural Heart,
Medtronic, BSCI

Support of this program is provided by Abbott, and the speaker is presenting on behalf of Abbott. This program is not intended for continuing education credits for any healthcare professional. All information presented is consistent with applicable FDA requirements.

All relevant financial relationships have been mitigated.

Faculty disclosure information can be found on the app

Navitor PAS Design

| | |
|----------------------------|--|
| Design | The Navitor Post-Approval Study (PAS) collected commercial use data through the TVT Registry |
| Clinical Sites | 2,958 patients in 198 centers within the United States |
| Population (N=2958) | Patients with tricuspid aortic valve morphology undergoing native TAVR between January 2023 - December 2024 |
| Key analyses | <ul style="list-style-type: none">• Present RWE with Navitor Classic and Navitor Vision• Technical success• Death or stroke at 30 days |

Study sponsored by Abbott.

The views or opinions presented here do not represent those of the American College of Cardiology Foundation, The Society of Thoracic Surgeons, or the STS/ACC TVT Registry. The industry data file upon which the analysis was performed is a subset of the full STS/ACC TVT Registry data submissions.



TCT

Garcia, Santiago. Real World Experience with the Navitor Valve in US Patients. New York Valves, June 2025, NY, USA.

See Important Safety Information referenced within.

MAT-2512512 v1.0 | Item approved for U.S. use only.

Baseline Characteristics

| Patient Characteristic | N=2958 | Medical History | N=2958 |
|--------------------------|-------------|---------------------------|--------|
| Age (years) | 81.4 ± 7.6 | A fib / flutter | 39.7% |
| Female | 62.4% | Conduction defect | 32.5% |
| BMI (kg/m ²) | 28.6 ± 6.9 | CABG | 11.0% |
| STS Score | 6.6 ± 5.4 | PCI | 29.4% |
| NYHA III or IV | 57.7% | Diabetes | 38.4% |
| Creatinine (mg/dL) | 1.3 ± 1.2 | Peripheral artery disease | 17.1% |
| LVEF (%) | 58.5 ± 11.0 | Permanent pacemaker | 12.1% |

Data presented as % or mean ± standard deviation.



TCT

Garcia, Santiago. Real World Experience with the Navitor Valve in US Patients. New York Valves, June 2025, NY, USA.

See Important Safety Information referenced within.

MAT-2512512 v1.0 | Item approved for U.S. use only.

Procedural Outcomes

| Outcomes | All Navitor (N=2958) | Navitor Classic (n=2237) | Navitor Vision (n=721) |
|-----------------------------|-------------------------|-----------------------------|---------------------------|
| Technical success* | 97.9% | 97.9% | 98.1% |
| Second valve required | 0.1% | 0.1% | 0.1% |
| Procedural mortality | 0.1% | 0.1% | 0.1% |
| Conversion to surgery | 0.1% | 0.1% | 0.0% |
| Major vascular complication | 1.7% | 1.8% | 1.4% |

*Technical success (upon leaving procedure room) defined as:

1. Freedom from mortality
2. Successful access, delivery of the device, and retrieval of the delivery system
3. Correct positioning of a single Navitor valve
4. Freedom from surgery or intervention related to the device or to a major vascular or access-related, or cardiac structural complication



TCT

Garcia, Santiago. Real World Experience with the Navitor Valve in US Patients. New York Valves, June 2025, NY, USA.

See Important Safety Information referenced within.

MAT-2512512 v1.0 | Item approved for U.S. use only.

30-Day Safety Events

| Outcomes | All Navitor (N=2958) | Navitor Classic (n=2237) | Navitor Vision (n=721) |
|---|-------------------------|-----------------------------|---------------------------|
| All-cause mortality or stroke | 5.2% | 5.2% | 4.9% |
| All-cause mortality | 2.8% | 2.9% | 2.5% |
| Stroke | 2.7% | 2.6% | 2.7% |
| Life-threatening or major bleeding | 0.7% | 0.8% | 0.6% |
| Major vascular complication | 1.8% | 1.9% | 1.5% |
| Aortic valve reintervention | 0.3% | 0.4% | 0.1% |
| New permanent pacemaker | 17.8% | 17.8% | 17.8% |
| Patients w/o baseline conduction defect | 14.1% | 14.5% | 12.3% |

Conduction defect defined as right or left BBB, sick sinus syndrome, or first-, second-, or third-degree heart block



TCT



Garcia, Santiago. Real World Experience with the Navitor Valve in US Patients. New York Valves, June 2025, NY, USA.

See Important Safety Information referenced within.

MAT-2512512 v1.0 | Item approved for U.S. use only.

PPI by Site Experience

Significant improvement in PPI rate after 20 implants



Number of procedures reflects procedures documented in TVT Registry from January 2023 to December 2024




TCT

Garcia, Santiago. Real World Experience with the Navitor Valve in US Patients. New York Valves, June 2025, NY, USA.
See Important Safety Information referenced within.

MAT-2512512 v1.0 | Item approved for U.S. use only.


Echocardiographic Outcomes

Hemodynamics



Site-assessed.

Paravalvular Leak



TCT


Garcia, Santiago. Real World Experience with the Navitor Valve in US Patients. New York Valves, June 2025, NY, USA.
See Important Safety Information referenced within.

MAT-2512512 v1.0 | Item approved for U.S. use only.

Four-Year Outcomes in the Navitor IDE Study

Device Durability

Bioprosthetic Valve Dysfunction



Component Rates at 4 Years

| | |
|--|-------------|
| Bioprosthetic Valve Dysfunction | 5.9% |
| Moderate HSVD | 0% |
| Non-structural valve deterioration | 1.7% |
| ➤ Severe PPM | 1.7% |
| ➤ Severe PVL | 0% |
| Infective endocarditis | 4.2% |
| Clinical valve thrombosis | 0% |
| Bioprosthetic Valve Failure | 0% |
| Severe HSVD | 0% |
| Aortic valve reintervention | 0% |
| Valve-related death | 0% |

Kaplan-Meier rates shown.



TCT

Manoharan, Ganesh. Four-Year Outcomes in the Navitor IDE Study. TCT, October 2025, San Francisco, CA, USA.

See Important Safety Information referenced within.

MAT-2512512 v1.0 | Item approved for U.S. use only.

Conclusions

- First report of commercial experience with Navitor Vision Valve
- Largest reported cohort of Navitor patients (N=2958)
- High technical success rate and low clinical event rates at 30 days
- Excellent valve hemodynamics and low rates of PVL
- Significant decrease in rates of PPI with increased device experience

Navitor and Navitor Vision offer a safe and effective treatment option for patients with severe symptomatic AS



TCT

Garcia, Santiago. Real World Experience with the Navitor Valve in US Patients. New York Valves, June 2025, NY, USA.
See Important Safety Information referenced within.

MAT-2512512 v1.0 | Item approved for U.S. use only.