

TCT 2025

***TAVIPILOT –
A unique AI&Robotic solution
for optimizing TAVI Procedures***

Mircea Moscu, Caranx Medical (CarvOlix Group), PhD




TCT®

TRANSCATHETER
CARDIOVASCULAR
THERAPEUTICS®


Disclosure of Relevant Financial Relationships

I, [Mircea Moscu](#), am presenting as Project Lead with Caranx Medical.

Why Robotics and AI?



Access to TAVI and outcomes can still be improved



1000'

*Thousands of patients eligible to TAVI remain untreated worldwide due to **lack of operators***

~10%

*Due to THV depth, **~10% of patients have conduction disorders and need pacemaker post TAVI***¹

~3%

*Due to THV depth, **~3% of patients have stroke post TAVI***¹

~2x

*TAVI centers performing **less than 100/year** have a **~2X higher mortality** rate than other TAVI centers¹*

75%

75% of cardiologists identify valve positioning as the most critical steps, followed by valve delivery²

TAVIPILOT Solution

TAVIPILOT Software – FDA cleared:
Real-time intra-operative guidance of TAVI
with millimetric precision


TAVIPILOT Robot – in development,
expecting FDA clearance 2026

TAVI Catheter Driver:

- Potentially reducing to **1 single operator**
- Driven by **TAVIPILOT Software**

TAVIPILOT solution aims to:

- **Increase accuracy** for valve positioning (mm precision)
- Reach **less variability** between users
- Potentially **reduce side-complications** such as pacemaker implant, etc.



TAVIPILOT Software is able to track anatomy and instruments in real-time on fluoroscopy images

- The AI detects and tracks the anatomy and instruments, following breathing and heart motion
- AI trained on the world's largest TAVI database (> 5000 patients)
- Once the contrast injected, the AI overlays NCC and initiates anatomical tracking
- Once the contrast dissipates, augmented reality continues tracking
- The implantation depth is measured for accurate positioning




TAVIPILOT Software: agnostic to C-arm imaging device

Siemens Artis

GE Discovery IGS7

Philips Azurion

Device



Compatibility



TAVIPILOT Robot (in development) – Principle of operation




Stent
positioning
phase




- **Foot pedal** in development could allow **single operator** use

TAVIPILOT Robot (in development)




- **Dedicated cassette** for each delivery device
- **Balloon or self-expandable valve** are managed manually by the operator
- **Compatible** with existing TAVI devices

Combined TAVIPILOT Software + Robot: Augmented teleoperation (in development)



- **AI** controls the robot
- **Clinician** controls the AI
- **Clinician** can revert at any time

AI Augmented Teleoperation: a faster learning curve and improved performance, tested on phantom



Towards autonomous robot-assisted transcatheter heart valve implantation: in vivo teleoperation and phantom validation of AI-guided positioning

Jonas Smits^{1*}, Pierre Schegg^{1†}, Loic Wauters¹, Luc Perard¹, Corentin Languepin¹, Davide Recchia¹, Vera Damerjian Pieters¹, Stéphane Lopez², Didier Tchetché³, Kendra Grubb⁴, Jorgen Hansen¹, Eric Sejor^{1,5} and Pierre Berthet-Rayne^{1,6}

TYPE Original Research
PUBLISHED 21 October 2025
DOI 10.3389/frobt.2025.1650228

- Tests performed by 3 TAVI experts
- 60 procedures on each test

Conclusion

- **TAVIPILOT Software: FDA cleared, the first real-time AI solution for intraoperative TAVI guidance** with millimetric precision
- **TAVIPILOT Robot:** In development, the first robotic solution for TAVI to facilitate valve positioning
- **TAVIPILOT Augmented Teleoperation:** Enhancing valve placement precision and democratizing TAVI

Please visit our booth #2252

Caranx Medical (CarvOlix group)



TCT®

TRANSCATHETER
CARDIOVASCULAR
THERAPEUTICS®