

# 2 Years Outcomes of Transfemoral J-VALVE for Chronic Aortic Regurgitation: A Prospective, Multicenter Study in 127 Cases

Jian'an Wang, MD

On Behalf of the J-VALVE TF China Investigators



**TCT**®


TRANSCATHETER  
CARDIOVASCULAR  
THERAPEUTICS®

# Disclosure of Relevant Financial Relationships

I, [Wang Jian'an, MD ] DO NOT have any financial relationships to disclose.


# Background

## Aortic Regurgitation




- No calcification
- Lack of anchoring area
- Annulus dilation

## Worse prognosis





# J-VALVE TF System



The valve sizes range from 21mm to 34 mm  
accommodating annular perimeters from 53 to 104 mm

# The Key J-VALVE TF Implantation Steps



## J-VALVE Positioning

Align with aortic annulus



## Anchor Deployment

Anchor ring into native cusps



## Valve Releasing

Auto-commissural alignment & open cell design benefit low coronaries

# **Trial Purpose**

**To evaluate the effectiveness and safety of the J-VALVE transfemoral aortic valve system amongst patients with symptomatic severe aortic regurgitation who are high-risk or inoperable for SAVR**

# Trial Design

**Prospective, Multicenter, Single Arm Evaluation of Patients with Symptomatic  $\geq 3+$  Aortic Regurgitation at High Risk or Inoperable for SAVR**

**J-VALVE TF Implantation**

**Clinical Evaluation, Echocardiography, NYHA and KCCQ etc  
at 30 Days, 6 Months, 1 Year and Annually up to 5 Years**

**30 Days Outcome**

**PCR London Valve 2024**

**1 Year Outcome**

**EuroPCR 2025**

**Comparison with Prespecified  
Performance Goal**

**2 Years Outcome**

**TCT 2025**

# Key Inclusion & Exclusion Criteria

## Inclusion:

- **Age  $\geq 65$  years;**
- Patients with symptomatic moderate to severe or severe aortic valve regurgitation, and NYHA  $\geq$  II
- High risk or inoperable for SAVR evaluated by the surgical team
- Aortic valve anatomy is suitable for TAVR evaluated by the investigators
- Sign informed consent form, and are willing to accept relevant examinations and clinical follow-ups



# Key Inclusion & Exclusion Criteria

## Exclusion:

- Acute myocardial infarction or coronary revascularization occurred within one month before procedure
- Cerebrovascular accident (CVA) occurred within 30 days before procedure
- Other valve diseases that need interventions;
- Previous aortic valve implantation (mechanical or biological)
- Left ventricular ejection fraction < 20%

# Primary Endpoints


**The primary endpoint was cumulative all-cause mortality at 12 months**

- All-cause mortality included cardiovascular mortality and non-cardiovascular mortality

# Key Secondary Endpoints

- Cardiovascular mortality
- Permanent pacemaker implantation
- Hemodynamic valve performance
- LV remodeling measured by echocardiography
- Functional Improvement of heart (NYHA)
- Quality of life (KCCQ)

# Screening and Patient Disposition



# Baseline Patient Characteristics

Variable	% or mean $\pm$ SD
Age (years)	73.9 $\pm$ 5.9
Female	36.2%
Mean STS Score	6.1 $\pm$ 4.5
NYHA Class III or IV	74.0%
Coronary artery Disease	45.7%
Frailty	74.0%
Hypertension	80.3%
Diabetes	11.8%

Variable	%
Prior Permanent Pacemaker	1.6%
Left Bundle Branch Block	7.1%
Right Bundle Branch Block	6.3%
Renal Insufficiency	12.6%
Pulmonary hypertension	15.7%
Peripheral Arterial Disease	58.3%
Atrial fibrillation	18.9%
Prior CVA or TIA	15.7%

# Baseline Echo Characteristics

Variable	% or mean $\pm$ SD
AR Severity	
Severe	<b>78.7%</b>
Moderate to Severe	21.3%
Pure AR	<b>89.0%</b>
AR with mild AS	11%
Vena Contracta Width (mm)	<b>7.5<math>\pm</math>1.7</b>
Mean Gradient (mmHg)	13.8 $\pm$ 5.0

Variable	% or mean $\pm$ SD
Ascending Aortic Diameter (mm)	40 $\pm$ 4.2
Mitral regurgitation (mild)	44.9%
Mitral regurgitation ( $\geq$ moderate)	20.5%
LVEDD (mm)	41.5 $\pm$ 8.8
LVEDD (mm)	<b>59.5<math>\pm</math>7.3</b>
LVEF	56.6 $\pm$ 11.3
PASP (mmHg)	32.8 $\pm$ 9.8

# Baseline CT Characteristics

Variable	% or mean $\pm$ SD
Leaflet	
Tricuspid	<b>96.1%</b>
Bicuspid/Quadricuspid	3.9%
Annular perimeter (mm)	<b>81.3<math>\pm</math>6.9</b>
> 80mm	<b>62.2%</b>
Leaflet or annular calcification	
No calcification	<b>76.4%</b>

Variable	% or mean $\pm$ SD
Leaflet or annular calcification	
Mild calcification	22.1%
LCA Height (mm)	12.8 $\pm$ 3.5
RCA Height (mm)	16.7 $\pm$ 3.9
Mean aortic annular angle ( $^{\circ}$ )	55.5 $\pm$ 10.9
> 70 $^{\circ}$ (%)	10.2%
Dextrocardia (%)	0.8%

# Procedural Outcomes

Outcome	%	Outcome	%
In-procedural Death	0%	Valve thrombosis	0%
Stroke	0%	Mitral valve damage or dysfunction	0%
Acute myocardial infarction	0%	Cardiac tamponade	0%
Bleeding	0%	Endocarditis	0%
Acute kidney injury	0%	Ventricular Perforation	0%
Converted to SAVR	2.4%	Aortic Dissection	0%
Valve in Valve	3.9%	Annular Rupture	0%
Coronary Obstruction	0%	Technical Success*	93.7%



# Safety Outcomes

Safety Outcome	30 Days	1 Year	2 Years
All cause mortality	1.6%	3.2%	6.3%
Cardiovascular mortality	1.6%	2.4%	3.9%
New permanent pacemaker implantation	9.5%	12.6%	13.4%
III ° AV block	3.9%	5.5%	5.5%
Major Vascular Complication	0.8%	1.6%	3.2%
Myocardial infarction	0%	0%	0%
All Stroke	0%	2.4%	5.5%
Major bleeding (life-threatening or disabling)	0%	0.8%	2.4%
Acute kidney injury	0%	0.8%	1.6%

# Cause of Death


Death 1Y	Day	CEC Adjudication
Aortic Dissection	11	Cardiovascular
Sudden Death	17	Cardiovascular
Hypertension, Heart Failure	139	Cardiovascular
Unknown	351	Non-Cardiovascular


Death 2Y	Day	CEC Adjudication
Aortic Dissection	391	Cardiovascular
Hemorrhagic stroke	423	Non-Cardiovascular
Sudden Death	468	Cardiovascular
Ischemic stroke	503	Non-Cardiovascular

# Hemodynamic Valve Performance



$p < 0.001$




# Paravalvular Regurgitation




# Left Ventricular Remodeling



# NYHA Improvement



# KCCQ Improvement



# Conclusion

**The Transfemoral J-VALVE has demonstrated the following characters in AR patients:**

- Low mortality and morbidity
- Low new permanent pacemaker implantation rate
- Excellent hemodynamics valve performance
- Echocardiography demonstrated significant improvement in LV remodeling
- Significant clinical functional improvement



# Acknowledgement

The study principal investigators and sponsor express our deepest gratitude to all the clinical research assistants, clinical research coordinators, CEC members, and all the patients at J-VALVE TF Trial in China