

# The Splitter: Enabler of transcatheter aortic valve-in-valve implantation

Yair Feld M.D., Ph.D., Arthur Kerner M.D

---



**TCT**®

TRANSCATHETER  
CARDIOVASCULAR  
THERAPEUTICS®

# Disclosure of Relevant Financial Relationships

Within the prior 24 months, I have had a financial relationship with a company producing, marketing, selling, re-selling, or distributing healthcare products used by or on patients:

## Nature of Financial Relationship

Consultant Fees/Honoraria

Individual Stock(s)/Stock Options

Royalties/Patent Beneficiary


## Ineligible Company

Medtronic, Edwards(s)


Cathalert, Cuspa, HVT, Paragate

HVT

# US ViV Market Forecast




***Redo-TAV replacement with S3-in-S3 and Evolut-in-S3 associated in 40% with intermediate or high risk for coronary obstruction***



***Miho Fukui. Circulation: Cardiovascular Interventions. Feasibility of Redo-Transcatheter Aortic Valve Replacement in Sapien Valves Based on In Vivo Computed Tomography Assessment, Volume: 16, Issue: 11, Pages: e013497, DOI: (10.1161/CIRCINTERVENTIONS.123.013497)***


# S3 in Evolut 20-70% high risk for coronary obstruction



***Feasibility of redo-TAVI in self-expanding Evolut valves: a CT analysis from the Evolut Low Risk Trial substudy" Euro-Intervention April 2023***






# Leaflet modification (BASILICA)




# Mechanical Vs. Electrosurgical

*ShortCut Catheter / Pi Cardia / UNICORN*



# Making a cut is not enough




- **Narrow Splay with newer devices**
- **Commissures need to be aligned**

**Khan et al. JACC Cardiovasc Interv. 2020 March 23; 13(6): 787–789.**




# Coronary obstruction is not the whole story

What about Coronary access?




# The Splitter

- The Splitter device performs cusp splitting and partial leaflet excision
- The device creates an intentional excision of valve's cusp tissue by using a steerable catheter with an electro-cutting wire loop running through a cutting head that mimics alligator jaws
- The Splitter is advanced over a standard 0.035" stiff guidewire, hence it is streamlined and integrated into the workflow of the valve implantation procedure.



# The Splitter - Animation



# Bench Test – Cutting Tissue

- Precise excision of leaflet tissue and creation of a large u-shaped window in the cusp
- The excised tissue is trapped inside the cutting head and removed through the catheter



The excised segment




The excised tissue trapped in the cutting head






The excised tissue removed

Human heart, ex-vivo









# Bench Test - Silicone Heart Model

The Splitter inside silicone model of human heart –  
access to valve plane and distal tip maneuvering



***Cutting head***



**Aortic arch**




**Cutting head**

**Control  
handle**

# Pre-Clinical Studies: Pig Study

- Device steering to the aortic valve
- Grasping the selected leaflet
- Laceration and removal of the leaflet segment
- Procedure done under fluoroscopy and echocardiography imaging guidance

*Distal end of the Splitter device*





*Cutting head*






# Pre-Clinical Studies: Pig Study – Fluoroscopy Imaging Movies






# R&D Animal Studies: Pig Study on an Implanted Prosthetic Valve

## Fluoroscopy imaging



**Cutting head  
exists the  
external unit**

## Cine imaging




# R&D Animal & Bench Studies: on a Prosthetic Valve

Prosthetic valve\* -  
post Splitter cutting  
(bench test)



Prosthetic valve\* -  
post Splitter cutting, in-vivo  
(explanted from pig post animal study)




\*Tested on Medtronic Evolut FX, 29 mm valve


HVT MEDICAL

# The Splitter – Essential to Remember

- Benefits:
- Precise control of the cutting site
- During the excision, the grasping jaws are stationary and stable (no pulling forces) while a running wire with electrical energy performs precise excision of a horseshoe shape
- Confirmation of a successful cutting
- Wide splay
- Electrosurgical cutting
- Fully integrated into the workflow of the TAVI procedure



# The Splitter – Status and Expected Timeline



## ➤ Splitter improvements:

### ➤ Improved steering

### ➤ Improved Positioning & Grasping

An adaptive jaw capable of dynamic grasping of both normal & calcified leaflets.

Verify the optimal leaflet position inside the Splitter jaws

### ➤ Improved cutting



# *Thank You*

*Yair Feld M.D., Ph.D.*

*Director – Structural and Interventional Cardiology*

*Tzafon Medical Center*

*Israel*

