

A Novel Technique for Commissural Alignment of Balloon-Expandable Transcatheter Aortic Valve Replacement: The SBCA Experience

Michael Shenoda, MD, FACC, FSCAI

Santa Barbara Cottage Hospital

Structural Heart Program



TCT®

TRANSCATHETER
CARDIOVASCULAR
THERAPEUTICS®

Disclosure of Relevant Financial Relationships

Within the prior 24 months, I have had a financial relationship with a company producing, marketing, selling, re-selling, or distributing healthcare products used by or on patients:

Nature of Financial Relationship

Consultant Fees/Honoraria

Ineligible Company

Edwards Lifesciences

Background

- Commissural alignment (CA) of balloon-expandable transcatheter aortic valve replacement (BE-TAVR) remains elusive
- Despite previous attempts of CA of BE-TAVR, none have been successful at consistently achieving CA
- Proposed benefits of CA include simpler coronary re-access, improved hemodynamics, decreased incidence of hypoattenuated leaflet thickening, and decelerated TAVR deterioration



Methods

- *Prospectively* evaluated the effects of BE-TAVR crimping orientation based on CT-derived aortic valve plane for *76 patients*, utilizing the **Shenoda-Blocker Commissural Alignment (SBCA)** technique
- Rates of successful CA were then *retrospectively* compared the previous *76 patients* who had undergone BE-TAVR using standard crimping and valve implantation techniques.

Methods

- Used the **ALIGN-TAVR Consortium** definitions of CA.
- However, valve alignment was assessed by *intraprocedural transesophageal echocardiogram* during the index TAVR, rather than retrospective CT analysis.
- Categories included:
 - CA (0-15.0°)
 - Mild commissural malalignment (CMA) (15.1°- 30.0°)
 - Moderate CMA (30.1°- 45.0°)
 - Severe CMA (45.1°- 60.0°)

CT-Derived Aortic Root Angle



SBCA Commissural Alignment Technique

CT-DERIVED AORTIC ROOT ANGLE	CRIMPING ORIENTATION
$< 35^\circ$	3:00
$35^\circ - 55^\circ$	3:30
$55^\circ - 65^\circ$	4:00
$> 65^\circ$	4:30

SBCA Commissural Alignment Technique




Results

- Baseline patient characteristics for both groups were similar
- Rate of CA for BE-TAVR utilizing the **SBCA technique** was **66%** (50/76 *patients*) compared to **31%** in the **standard group** (23/75 *patients*)
- *Two-Sample Proportions Chi-Squared Test* was performed to examine the relationship between these two independent proportions, [$\chi^2(1, N=151)=17.27, 95\% \text{ CI } (0.18, 0.51), p < 0.01$] signifying statistically sufficient evidence to conclude that the proportion of patients with CA utilizing the SBCA technique is higher than the proportion of patients with CA pre-technique

Results

- The rate of **CA + mild CMA (coronary alignment)** for the SBCA technique was **80%** (61/76 patients) compared to **53%** (40/75 patients) in the **standard group**, which was also statistically significant for the SBCA technique (95% CI (0.11, 0.43), $p < 0.01$).
- **Severe CMA** was rare with the SBCA technique at **1.3%** (1/76 patients) compared to the **standard technique** at **31%** (23/75 patients) favoring the SBCA technique (95% CI (0.17, 0.41), $p < 0.01$).



Aortic Valve Plane Trends



Bicuspid Valves


- Bicuspid type 0: 2 (1 Mod CMA, 1 Severe CMA)
- Bicuspid type 1: 4 (2 CA, 1 Mild CMA, 1 Mod CMA)

Example – Trileaflet Valve




37° Aortic Valve Plane

Example – Trileaflet Valve



37° Aortic valve plane – 3:30 crimping
orientation

Example – Bicuspid Valve Type I





57° Aortic Valve Plane

Example – Bicuspid Valve Type I




57° Aortic valve plane - 4:00 crimping
orientation

Caution/Downfalls – Iliofemoral Tortuosity



Caution/Downfalls – Aortic Tortuosity



Caution/Downfalls – Catheter Rotation



- Balloon port orientation should always remain in the 9:00 position
- May need to adjust after crossing the valve
- Adjusting during rapid ventricular pacing results in better translation of torque

Conclusion

- SBCA technique of Sapien 3 THV orientation during valve crimping based on CT-derived aortic valve plane resulted in a significant improvement in **commissural alignment (66%)** when compared to **standard implantation technique (31%)**
- Increased to **80%** with **CA + mild CMA (coronary alignment)** compared to **53%** with the **standard implantation technique**
- Severe CMA was rare with the SBCA technique
- CA outliers were associated with *extremes of aortic root angulation ($<35^\circ$ or $>65^\circ$), significantly tortuous iliofemoral vessels, aortic tortuosity*, and with significant catheter manipulation during TAVR.

Citations

- Tang, GHL et al. Alignment of transcatheter aortic-valve neo-commissures (ALIGN TAVR): Impact on final orientation and coronary artery overlap. *JACC Cardiovasc Inter.* 2020, 13(9):1030-42
- Tang, GHL et al. Rationale, definitions, techniques, and outcomes of commissural alignment in TAVR; From the ALIGN-TAVR consortium. *JACC Cardiovasc Inter.* 2022, 15(15):1497-1518
- Akodad M et al. Transcatheter heart valve commissural alignment: An updated review. *Front. Cardiovasc. Med.* 2023, 10:1154556. doi: 10.3389/fcvm.2023.1154556