

# Optimizing TAVI: Utilizing Standardized Invasive Hemodynamics

*Ravi Ramana, DO, FACC*

*Interventional Cardiologist*

*Heart Care Centers of Illinois*



**TCT**<sup>®</sup>

TRANSCATHETER  
CARDIOVASCULAR  
THERAPEUTICS<sup>®</sup>

# Disclaimer

The speaker's presentation today is on behalf of Haemonetics Corporation. Any discussion regarding Haemonetics' products during the presentation today is limited to information that is consistent with the Haemonetics' labeling for those products. Please consult product labels and inserts for any indications, contraindications, hazards, warnings, precautions and instructions for use.

The clinical experiences, techniques, procedures and other medical/pharmaceutical product suggestions presented herein that are apart from Haemonetics' products are solely those of the presenter, not Haemonetics Corporation, and are offered for information purposes only. The results from these case studies/examples may not be predictive for all patients. Individual results may vary depending on a variety of patient-specific attributes. The physician/s have been compensated by Haemonetics Corporation for their time.

This presentation is intended for healthcare professionals only.


# >20 Years of TAVI: Advancing AS Management

Advances in TAVI & clinical evidence have transformed the treatment landscape and TAVR has become the standard of care




# >20 Years of TAVI: A Shift in Focus

- Advances in valve design, improved imaging and technical enhancements have contributed to increased safety and a decline in procedural complications
- The focus of TAVI is now shifting towards procedural efficiency and lifetime patient management



# >20 Years of TAVI: A Shift in Focus

- Advances in valve design, improved imaging and technical enhancements have contributed to increased safety and a decline in procedural complications
- The focus of TAVI is now shifting towards procedural efficiency and lifetime patient management



# But... are we as good ...

NEWS Conference News EACTS 2025  
**SAVR Bests TAVI for 5-Year Survival in Combined  
RCT, Observational Data**




But... are we as good as we think we are???



# Hemodynamic Performance

## *The importance of hemodynamics*

- *Observed discrepancies between echo and invasive gradients*
- *Invasive measurements are accurate but:*
  - *Take time to setup*
  - *Need wire and pigtail exchanges*
- *LV pressure is not monitored constantly during TAVI*
- *Large-amplitude LV pressure creates artifacts in fluid-filled lines*
- *Assessing and treating lower risk and younger patients*
- *Assessing and treating ASYMPTOMATIC AS*






# Hemodynamic Performance

## *The importance of hemodynamics*

- *Observed discrepancies between echo and invasive gradients*
- *Invasive measurements are accurate but:*
  - *Take time to setup*
  - *Need wire and pigtail exchanges*
- *LV pressure is not monitored constantly during TAVI*
- *Large-amplitude LV pressure creates artifacts in fluid-filled lines*
- *Assessing and treating lower risk and younger patients*
- *Assessing and treating ASYMPTOMATIC AS*




***We must not only be good, but we need to be nearly “perfect”***

# Echo Versus Invasive Gradients


## Echocardiographic Versus Invasive Aortic Valve Gradients in Different Clinical Scenarios

- *Echo is non-invasive but can be inaccurate*
- *Pre-procedure gradients/ViV-TAVI impact both echo/invasive post-procedure gradients*
- *TAVI valve types impact echo/invasive mean gradients differently*



Abbas, Amr E. et al. Echocardiographic Versus Invasive Aortic Valve Gradients in Different Clinical Scenarios. Journal of the American Society of Echocardiography, Volume 36, Issue 12, 1302 - 1314

# Echo-derived Gradients Do Not Correlate with Invasive Normal Functioning Transcatheter Heart Valves



- Post-TAVR echo overestimates transvalvular mean gradients compared with invasive measurements

And poor correlation suggests these modalities cannot be used interchangeably

# Echo-Derived Gradients Overestimate True Cath Gradient by ~1.6x for Balloon Expandable TAVR Valves




Table 1. Comparison of invasive transvalvular gradients of similar-sized balloon-expandable and self-expanding transcatheter heart valves immediately post implantation.

	N	CT Area (mm <sup>2</sup> )	Area-Derived Diameter (mm)	CT Perimeter (mm)	Perimeter-Derived Diameter (mm)	Invasive Gradient	Doppler Gradient	P-Value
<b>Small Valve</b>								
Total	128	370.7 ± 33.8	21.7 ± 1.0	66.6 ± 3.2	22.2 ± 1.0	7.7 ± 3.8	10.3 ± 5.4	<.001
BEV	61	371.8 ± 30.6	21.7 ± 0.9	69.8 ± 3.0	22.2 ± 0.9	8.1 ± 3.6	13.3 ± 5.3	<.001
SEV	67	369.7 ± 36.6	21.7 ± 1.1	69.5 ± 3.5	22.1 ± 1.1	7.4 ± 3.9	7.5 ± 4.0	.87
P-value		.74	.70	.57	.57	.31	<.001	
<b>Medium Valve</b>								
Total	240	464.8 ± 39.1	24.3 ± 1.0	77.9 ± 3.2	24.8 ± 1.0	7.9 ± 4.1	10.3 ± 4.1	<.001
BEV	164	467.1 ± 36.4	24.4 ± 0.9	78.1 ± 3.1	24.9 ± 1.0	7.6 ± 4.0	11.0 ± 3.8	<.001
SEV	76	459.9 ± 44.3	24.2 ± 1.2	77.4 ± 3.6	24.7 ± 1.2	8.7 ± 4.3	8.6 ± 4.2	.94
P-value		.26	.24	.22	.22	.05	<.001	
<b>Large Valve</b>								
Total	131	576.3 ± 49.3	27.1 ± 1.1	86.7 ± 3.6	27.6 ± 1.1	7.2 ± 3.0	9.2 ± 3.9	<.001
BEV	95	578.8 ± 48.7	27.1 ± 1.1	86.8 ± 3.6	27.6 ± 1.1	6.9 ± 3.0	10.0 ± 3.8	<.01
SEV	36	569.4 ± 51.1	26.9 ± 1.2	86.6 ± 3.5	27.6 ± 1.1	8.1 ± 2.8	7.2 ± 3.3	.21
P-value		.38	.37	.74	.74	.05	<.001	

BEV = balloon-expandable valve; CT = computed tomography; SEV = self-expanding valve.

# Echo-Derived Gradients Overestimate True Cath Gradient by ~1.6x for Balloon Expandable TAVR Valves

## No Difference in Invasive Gradients Between BEV and SEV




Table 1. Comparison of invasive transvalvular gradients of similar-sized balloon-expandable and self-expanding transcatheter heart valves immediately post implantation.

	N	CT Area (mm <sup>2</sup> )	Area-Derived Diameter (mm)	CT Perimeter (mm)	Perimeter-Derived Diameter (mm)	Invasive Gradient	Doppler Gradient	P-Value
<b>Small Valve</b>								
Total	128	370.7 ± 33.8	21.7 ± 1.0	66.6 ± 3.2	22.2 ± 1.0	7.7 ± 3.8	10.3 ± 5.4	<.001
BEV	61	371.8 ± 30.6	21.7 ± 0.9	69.8 ± 3.0	22.2 ± 0.9	8.1 ± 3.6	13.3 ± 5.3	<.001
SEV	67	369.7 ± 36.6	21.7 ± 1.1	69.5 ± 3.5	22.1 ± 1.1	7.4 ± 3.9	7.5 ± 4.0	.87
P-value		.74	.70	.57	.57	.31	<.001	
<b>Medium Valve</b>								
Total	240	464.8 ± 39.1	24.3 ± 1.0	77.9 ± 3.2	24.8 ± 1.0	7.9 ± 4.1	10.3 ± 4.1	<.001
BEV	164	467.1 ± 36.4	24.4 ± 0.9	78.1 ± 3.1	24.9 ± 1.0	7.6 ± 4.0	11.0 ± 3.8	<.001
SEV	76	459.9 ± 44.3	24.2 ± 1.2	77.4 ± 3.6	24.7 ± 1.2	8.7 ± 4.3	8.6 ± 4.2	.94
P-value		.26	.24	.22	.22	.05	<.001	
<b>Large Valve</b>								
Total	131	576.3 ± 49.3	27.1 ± 1.1	86.7 ± 3.6	27.6 ± 1.1	7.2 ± 3.0	9.2 ± 3.9	<.001
BEV	95	578.8 ± 48.7	27.1 ± 1.1	86.8 ± 3.6	27.6 ± 1.1	6.9 ± 3.0	10.0 ± 3.8	<.01
SEV	36	569.4 ± 51.1	26.9 ± 1.2	86.6 ± 3.5	27.6 ± 1.1	8.1 ± 2.8	7.2 ± 3.3	.21
P-value		.38	.37	.74	.74	.05	<.001	

BEV = balloon-expandable valve; CT = computed tomography; SEV = self-expanding valve.

# Why?:

## 2D Echo Methods of AVA Calculation may not account for Valve Design



$$AVA = \frac{\pi(\text{radius LVOT})^2 \times (V_{\max\text{LVOT}})}{(V_{\max\text{AV}})}$$

Actual Valve EOA

EOA Reported by ECHO

- Also, ECHO VELOCITY is a proxy for PRESSURE

1344

JACC Vol. 12, No. 5  
November 1988:1344-53

### Review of Hydrodynamic Principles for the Cardiologist: Applications to the Study of Blood Flow and Jets by Imaging Techniques

AJIT P. YOGANATHAN, PhD, EDWARD G. CAPE, BS, HSING-WEN SUNG, PhD,  
FRANK P. WILLIAMS, PhD, ABDUL JIMOH, MS  
Atlanta, Georgia


**When does the simplified Bernoulli equation “not work”?**  
(i.e.,  $\Delta P_{\text{peak}} = 4V_2^2$ ). 1) When the proximal velocity is of the same order of magnitude as the distal velocity (1–3). Examples are: a) aortic regurgitation in combination with aortic stenosis; and b) prosthetic heart valves. Note, that a 1 to 2 m/s proximal velocity leads to a 4 to 16 mm Hg decrease in pressure gradient. In such cases use equation 3, that is, the Bernoulli equation.

# Why?:

## Flow Pattern Influence how Echo Measures Velocity

### Laminar flow

An **efficient** design results in laminar flow and maintained velocity between the LV and the aorta.




Laminar flow =  
higher velocity

**Higher** echo  
derived  
gradient



### Turbulent flow

An **inefficient** design results in turbulent flow and inconsistent velocity between LV and the aorta.



Turbulent flow =  
lower velocity

**Lower** echo  
derived  
gradient



Because laminar flow maintains velocity, **echo overestimates the pressure gradient in a valve with laminar flow** compared to a valve with turbulent flow perversely rewarding valve design with turbulent flow.

# What does all mean?

- I don't know....





# What does all mean?

- I don't know....



But, more data is a good thing.

# The Importance of Hemodynamics

## Prognostic value of invasive versus echocardiography-derived aortic gradient in patients undergoing TAVI

- 806 TAVI patients: EDW S3 Balloon Expandable Valves (BEV, n=474) or MDT Evolut Self-Expanding Valves (SEV, n=332) from 2014 to 2023
- Examined the prognostic value of invasive and echo-derived gradients after implantation of SEV and BEV and to compare gradients for SEV versus BEV

### Prognostic value of invasive versus echocardiography-derived aortic gradient in patients undergoing TAVI

Mark M.P. van den Dorpel<sup>1</sup>, MD; Stram Chatterjee<sup>1</sup>, BSc; Rik Adrichem<sup>1</sup>, MD; Sarah Verhemeij<sup>1</sup>, MD; Jaskella Kandy<sup>1</sup>, MD, PhD; Ranga-Jung Nui<sup>1</sup>, MD, PhD; Joost Daemen<sup>1</sup>, MD, PhD; Claire Ben Ben<sup>1</sup>, MD, PhD; Alexander Hirsch<sup>2</sup>, MD, PhD; Marcel L. Geleijne<sup>1</sup>, MD, PhD; Nicolas M. Van Mieghem<sup>1,\*</sup>, MD, PhD

<sup>1</sup>Corresponding author: Department of Cardiology, Thoraxcenter, Erasmus University Medical Center, office Nt-645, Dr. Molewaterplein 40, 3015 GD, Rotterdam, the Netherlands. E-mail: [n.vanmieghem@erasmusmc.nl](mailto:n.vanmieghem@erasmusmc.nl)

This paper also includes supplementary data published online at: <http://eurointervention.pcronline.com/doi/10.4244/EIJ-D-24-00341>

**BACKGROUND:** Recent studies report a discordance between invasive and echocardiography-derived gradients after transcatheter aortic valve implantation (TAVI) with balloon-expandable (BEV) and self-expanding valves (SEV). There are limited data on the determinants and clinical implications of this discordance.

**AIMS:** We aimed to examine the prognostic value of invasive and echocardiography-derived gradients after implantation of SEV and BEV and to compare gradients for SEV versus BEV.

**METHODS:** We performed a retrospective, propensity score-matched study. Invasive measurements were obtained before and immediately after TAVI. Echocardiography was performed before and within 24 hours after TAVI, and at 1 year. Clinical outcomes were assessed at 30 days, 1 year, and 2 years.

**RESULTS:** The 1:1 propensity score matching resulted in 436 matched pairs (436 SAPIEN 3 and 436 Evolut). Invasive gradients post-TAVI independently predicted higher risk for all-cause mortality at 30 days, 1 year and 2 years as a continuous variable (hazard ratio [HR] 1.07, 95% confidence interval [CI]: 1.00-1.14; p=0.038; HR 1.06, 95% CI: 1.01-1.11; p=0.007; HR 1.05, 95% CI: 1.01-1.09; p=0.011, respectively) and by using >10 mmHg as a cutoff (HR 1.95, 95% CI: 1.13-3.78; p=0.028; HR 1.91, 95% CI: 1.11-3.65; p=0.030; HR 1.61, 95% CI: 1.03-2.96; p=0.021, respectively), but echocardiography-derived gradients did not (HR 1.13, 95% CI: 0.87-1.75; p=0.247; HR 1.02, 95% CI: 0.95-1.10; p=0.639; HR 0.99, 95% CI: 0.94-1.07; p=0.979, respectively). Mean gradients before and after TAVI were higher by echocardiography than by invasive measurements. The difference was more pronounced after implantation with BEV than SEV (7.0 [25th-75th percentile: 4.0-11.0] mmHg vs 5.0 [2.0-7.0] mmHg; p=0.001). Smaller valve size, higher ejection fraction and higher stroke volume amplified the discordance. Invasive mean gradients were similar after SEV and BEV (3.0 [0.0-6.0] mmHg vs 3.0 [0.0-6.0] mmHg; p=0.166), but echo-derived mean gradients were lower after SEV versus BEV (6.0 [6.0-11.0] mmHg vs 11.0 [8.0-14.0] mmHg; p=0.001).

**CONCLUSIONS:** Only invasively measured but not echocardiography-derived transvalvular mean gradients correlate with 30-day, 1-year and 2-year mortality. Aortic gradient measurements are higher by echocardiography than by invasive assessment and more so for BEV than SEV. Smaller valve size, higher ejection fraction and higher stroke volume increase this discordance between echocardiography and invasive assessment.

**KEYWORDS:** balloon-expandable valve; echocardiography; invasive; pressure gradient; self-expanding valve; TAVI

# The Importance of Hemodynamics


## Prognostic value of invasive versus echocardiography-derived aortic gradient in patients undergoing TAVI

### Echo overestimates gradients

- **SEV:** Echo pressure gradient (PG) 8.0 mmHg vs Invasive PG 3.0 mmHg
- **BEV:** Echo PG 11.0 mmHg vs Invasive PG 3.0 mmHg

### Contributing factors

- Small valve size
- High LVEF
- High stroke volume



# Independent Predictors of Mortality

**Table 5. Cox proportional hazards regression model for all-cause mortality.**

Univariable	30 days		1 year		2 years	
	HR (95% CI)	p-value	HR (95% CI)	p-value	HR (95% CI)	p-value
<b>Echo gradients</b>						
Echo mean gradient post-TAVI (continuous)	1.13 (0.87-1.75)	0.248	1.02 (0.95-1.10)	0.639	0.99 (0.94-1.07)	0.979
Echo mean gradient post-TAVI >10 mmHg	1.07 (0.95-1.19)	0.189	1.76 (0.51-7.66)	0.309	1.32 (0.66-11.02)	0.573
Echo mean gradient post-TAVI >20 mmHg	1.05 (0.95-1.14)	0.205	2.37 (0.33-17.23)	0.394	2.33 (0.32-16.90)	0.404
<b>Invasive gradients</b>						
Invasive mean gradient post-TAVI (continuous)	1.08 (1.01-1.15)	0.034 <sup>&amp;</sup>	1.05 (1.00-1.10)	0.047 <sup>&amp;</sup>	1.05 (1.01-1.09)	0.009 <sup>&amp;</sup>
Invasive mean gradient post-TAVI >10 mmHg	2.94 (1.12-8.31)	0.030 <sup>&amp;</sup>	1.83 (1.07-3.31)	0.042 <sup>&amp;</sup>	1.71 (1.04-2.82)	0.036 <sup>&amp;</sup>

# Independent Predictors of Mortality

**Table 5. Cox proportional hazards regression model for all-cause mortality.**

Univariable	30 days		1 year		2 years	
	HR (95% CI)	p-value	HR (95% CI)	p-value	HR (95% CI)	p-value
<b>Echo gradients</b>						
Echo mean gradient post-TAVI (continuous)	1.13 (0.87-1.75)	0.248	1.02 (0.51-2.10)	0.639	0.99 (0.94-1.07)	0.979
Echo mean gradient post-TAVI >10 mmHg	1.07 (0.95-1.19)	0.189	1.76 (0.51-7.66)	0.399	1.32 (0.66-11.02)	0.573
Echo mean gradient post-TAVI >20 mmHg	1.05 (0.96-1.14)	0.205	2.37 (0.33-17.23)	0.394	1.23 (0.32-16.1)	0.404
<b>Invasive gradients</b>						
Invasive mean gradient post-TAVI (continuous)	1.08 (1.01-1.15)	0.034 <sup>&amp;</sup>	1.07 (1.00-1.10)	0.047 <sup>&amp;</sup>	1.05 (1.01-1.09)	0.009 <sup>&amp;</sup>
Invasive mean gradient post-TAVI >10 mmHg	2.94 (1.12-8.31)	0.030 <sup>&amp;</sup>	1.83 (1.07-3.31)	0.042 <sup>&amp;</sup>	1.71 (1.04-2.82)	0.036 <sup>&amp;</sup>

# Independent Predictors of Mortality

**Table 5. Cox proportional hazards regression model for all-cause mortality.**

Univariable	30 days		1 year		2 years	
	HR (95% CI)	p-value	HR (95% CI)	p-value	HR (95% CI)	p-value
<b>Echo gradients</b>						
Echo mean gradient post-TAVI (continuous)	1.13 (0.87-1.75)	0.248	1.02 (0.95-1.10)	0.639	0.99 (0.94-1.07)	0.979
Echo mean gradient post-TAVI >10 mmHg	1.07 (0.95-1.19)	0.189	1.76 (0.51-7.66)	0.309	1.32 (0.66-11.02)	0.573
Echo mean gradient post-TAVI >20 mmHg	1.05 (0.96-1.14)	0.205	2.37 (0.51-11.0)	0.284	2.33 (0.51-10.9)	0.404
<b>Invasive gradients</b>						
Invasive mean gradient post-TAVI (continuous)	1.08 (1.01-1.15)	0.034 <sup>&amp;</sup>	1.05 (1.00-1.10)	0.047 <sup>&amp;</sup>	1.05 (1.01-1.09)	0.009 <sup>&amp;</sup>
Invasive mean gradient post-TAVI >10 mmHg	2.94 (1.12-8.31)	0.030 <sup>&amp;</sup>	1.83 (1.07-3.31)	0.042 <sup>&amp;</sup>	1.71 (1.04-2.82)	0.036 <sup>&amp;</sup>

# The Importance of Hemodynamics

## Conclusions:

### Prognostic value of invasive versus echocardiography-derived aortic gradient in patients undergoing TAVI

- **Invasive mean pressure gradients correlated with all-cause mortality**
  - At 30 days: HR 1.95 (p=0.028);
  - At 1 year: HR 1.91 (p=0.030);
  - At 2 years: HR 1.61

# Aortic Regurgitation Index - ARi

Journal of the American College of Cardiology  
© 2013 by the American College of Cardiology Foundation  
All rights reserved.

Vol. 61, No. 11, 2013  
ISSN 0885-0666  
DOI: 10.1016/j.jacc.2013.04.001

## CLINICAL RESEARCH


## Valvular Heart Disease

### Aortic Regurgitation Index Defines Severity of Peri-Prosthetic Regurgitation and Predicts Outcome in Patients After Transcatheter Aortic Valve Implantation

Jan-Malte Sinning, MD, Christoph Hammerstingl, MD, Mariuca Vasa-Nicotera, MD, Viktoria Adenauer, MD, Sisa Josefina Lema Cachiguango, MD, Anne-Catherine Scheer, MD, Sven Hausen, MD, Alexander Sedaghat, MD, Alexander Ghanem, MD, Cornelius Müller, MD, Eberhard Grube, MD, Georg Nickenig, MD, Nikos Werner, MD

Bonn, Germany

- 141 TAVR patients by ECHO, Angiogram and ARi
- $ARi = [(DBP - LVEDP) / SBP] \times 100$



**Figure 1** Calculation of the AR Index

Simultaneous determination of left ventricular end-diastolic pressure (LVEDP) (blue line) and diastolic blood pressure (DBP) in the aorta (red line) in a patient without peri-prosthetic aortic regurgitation (periAR) (A) and in a patient with moderate periAR (B) for the calculation of the aortic regurgitation (AR) index:  $[(DBP - LVEDP) / SBP] \times 100$ . (A) AR index =  $[(65 - 10) / 160] \times 100 = 34.4$ . (B) AR index =  $[(40 - 20) / 130] \times 100 = 15.4$ .




# Aortic Regurgitation Index - ARi

- In Aortic Regurgitation, blood flow backwards into the left ventricle: aortic diastolic pressure decreases while LVEDP increases


- ARi ↓ as AR/PVL ↑

## Balloon post-dilation for reduction of



$$ARi = 100 \times \frac{Ao \text{ diastolic} - LVEDP}{Ao \text{ systolic}}$$

# ARi Identifies Patients at Highest Risk of Mortality After TAVI



**Figure 4** 1-Year Mortality According to Severity of PeriAR and the AR Index


1-year all-cause mortality (%) according to the severity of periAR (none/mild vs. moderate/severe) and the AR index cutoff value. Abbreviations as in Figures 1 and 2.

- ***ARi > 25 → 1-yr mortality 15.4%***
- ***ARi ≤ 25 → 1-yr mortality 33.3% (2x higher)***
- ***Moderate/Severe periAR + ARi ≤ 25 → 1-yr mortality 70.6%***

# Aortic Regurgitation Index Ratio – ARi Ratio

- The ARi ratio, which reflects acute hemodynamic changes after TAVI, is useful to identify patients with negative outcomes<sup>1</sup>
- So, associated with  $\geq$  moderate PVL and increased one-year mortality:<sup>2</sup>
  - Post procedural ARi values  $<25$  AND
  - ARi ratio of  $<0.6$  (or  $<0.5$ )


## Calculation of the Aortic Regurgitation Index Ratio<sup>1</sup>



**Note:**  
Haemodynamic parameters such as ARi alone cannot provide basis for clinical decision and their interpretation is left to physicians. Values should be interpreted by comparing pre- to post procedural measurements, along with other information available such as imaging modalities.


$$\text{ARi ratio} = \frac{\text{ARi post}}{\text{ARi pre}}$$

# SavvyWire® Guidewire – Hemodynamics




## Unipolar left ventricular pacing, at anytime without catheter

- Built-in shaft insulation supports consistent capture and low thresholds
- Eliminates RV access for eligible patients
- Simply attach alligator clips to an external pulse generator



# SavvyWire® Guidewire – Hemodynamics


## Continuous measurement and display



- *Pulse rate*
- *Aortic pressures (aortic pigtail/transducer)*
  - ✓ *Systolic, diastolic*
- *Left ventricular pressures*
  - ✓ *Systolic, diastolic, LVEDP*
- *Transvalvular gradients*
  - ✓ *Mean, peak-to-peak, instantaneous*
- *Aortic Regurgitation indices*
  - ✓ *ARi, ARi ratio*

# SavvyWire® Guidewire – Hemodynamics


- Hemodynamics are always available, without exchange
  - At baseline
  - After pre-dilation
  - During partial deployment (self-expanding valves)
  - Before post-dilation
  - Final hemodynamic assessment
- Accurate, live diagnostics
  - Regurgitation – from diastolic gradient ( $Ao - LVEDP$ )



# Case Presentation


- 85yo F with Diastolic CHF, HTN, Frail reporting exertional dyspnea and edema. Progressive easy fatigue.
- ECHO revealed EF 65% with severe aortic stenosis (mAVG 38mmHg, Jet 4.2m/s, AVA 0.8cm<sup>2</sup>).
- Cath with mild non-obstructive CAD.
- CTA TAVR protocol suggest acceptable iliofemoral vascular.

# SavvyWire® Guidewire, During and Post Deployment






# SavvyWire® Guidewire, During and Post Deployment





# SavvyWire® Guidewire, Post-BAV Gradient and ARi



Post Dilate

# SavvyWire® Guidewire, Post-BAV Gradient and ARi



# SavvyWire® Guidewire – Hemodynamics

- What does it do for me?




# SavvyWire® Guidewire – Hemodynamics

- What does it do for me?
  - Continuous, accurate measurement of hemodynamic pressure
  - Delivers “gold standard” invasive pressure measurements
  - Minimized exchanges for pressure readings
  - “Immediate” hemodynamic assessment of residual AI/PVL
  - Evidenced-based, reliable pacing (one less vessel sheath)
  - Increased lab efficiency



# SavvyWire<sup>®</sup> Guidewire – Hemodynamics

- What does it do for me?
  - Continuous, accurate measurement of hemodynamic pressure
  - Delivers “gold standard” invasive pressure measurements
  - Minimized exchanges for pressure readings
  - “Immediate” hemodynamic assessment of residual AI/PVL
  - Evidenced-based, reliable pacing (one less vessel sheath)
  - Increased lab efficiency



Meaning: simpler, safer and accurate data

## *Science advances...*

Realistic expectations for life are that we are going to be better today than we were yesterday, be better tomorrow than we were today. That's a plan for success. So [the key is] simple: just work.

