

# RedoTAVR coronary obstruction risk in small annuli

*A post-TAVR CT study*

Gaetano Liccardo, MD. *ICPS, Massy, France*



**TCT**<sup>®</sup>

TRANSCATHETER  
CARDIOVASCULAR  
THERAPEUTICS<sup>®</sup>

# Disclosure of Relevant Financial Relationships

I, Gaetano Liccardo, DO NOT have any financial relationships to disclose.

# Background

- TAVR is a well-established treatment for severe AS.
- *Small aortic annuli* -> risk of PPM.
- SEVs in *small annuli*:
  - Superior hemodynamic performance/ less PPM.
  - Similar clinical outcomes.

TAVI is recommended in patients  $\geq 70$  years of age with tricuspid AV stenosis, if the anatomy is suitable.<sup>d 1-4,389-397,465,485,486</sup>

I

A

- The need for redoTAVR procedures is expected to grow in the coming years.

# Background - CO risk mechanism

An important consideration when planning a *redoTAVR* procedure is that the leaflets of the first THV will be displaced by the second prosthesis, forming a **neoskirt-covered stent**.




# Aim

To evaluate **CO risk** in case of **redoTAVR** in patients with *small* vs. *non-small* aortic annuli.



# Methods-CT analysis


- 167 post-TAVR CT scans were analyzed. Patients were stratified into small annuli ( $\leq 430 \text{ mm}^2$ ) and non-small annuli groups ( $> 430 \text{ mm}^2$ ).
- VTC and VTA were measured.



# Methods-CT analysis

- RedoTAVR CO risk evaluated at different risk planes. For SEVs, these planes corresponded to nodes 4, 5, and 6 on the valve frame. For BEVs, risk plane was placed at the outflow level of the device.
- High risk of CO when  $VTC < 4$  mm or  $VTA < 2$  mm (below the risk plane).







# THVs across the two groups

THV	Small Annuli (n, %)	Non-Small Annuli (n, %)
Sapien 3 Ultra 23 mm	42 (58.3%)	3 (3.2%)
Sapien 3 Ultra 26 mm	5 (6.9%)	49 (51.6%)
Sapien 3 29 mm	-	19 (20.0%)
Evolut Pro Plus 23 mm	3 (4.2%)	-
Evolut R/Pro Plus 26 mm	15 (20.8%)	1 (1.1%)
Evolut R/ Pro Plus 29 mm	7 (9.8%)	12 (12.6%)
Evolut Pro Plus 34 mm	-	11 (11.5%)

There was no significant association between annulus size and SEV use ( $\chi^2[1, N=167]=1.77$ ;  $p=0.184$ ).

# CO risk evaluation

Predicted CO risk at redoTAVR: overall and by annulus size

88/167 patients (53%) of the overall study population were considered at high of CO risk in case of redoTAVR, with no differences *between small or non-small annuli* (OR = 1.65, 95% CI: 0.89–3.06,  $p = 0.112$ ).

# CO risk: SEVs Node 6 vs BEVs

Compared to patients having received BEVs, those treated with a SEV had a higher post-redoTAVR risk of CO if they had a small (OR = 15.52, 95 % CI: 3.28-73.6,  $p < 0.001$ ), but not if they had a non-small (OR = 1.44, 95 % CI: 0.57-3.65,  $p = 0.441$ ) annulus.


# CO risk: SEVs Node 5 vs BEVs

Compared to patients having received BEVs, those treated with a SEV had a higher post-redoTAVR risk of CO if they had a small (OR = 3.13, 95 % CI: 1.13-8.71,  $p = 0.03$ ), but not if they had a non-small (OR = 0.73, 95 % CI: 0.28-1.89,  $p = 0.52$ ) annulus.

# CO risk: SEVs Node 4 vs BEVs

When the SEV risk plane was at node 4, SEVs were not associated with a higher post-redoTAVR risk of CO compared to BEVs in both patients with small (OR = 1.26, 95 % CI: 0.51-3.61,  $p = 0.54$ ) or non-small (OR = 0.73, 95 % CI: 0.28-1.88,  $p = 0.52$ ) annuli.

# CO risk: SEVs vs BEVs across the two groups



# Take Home Messages

- As the population undergoing TAVR has expanded to **younger and lower-risk** patients, the need for **redoTAVR** procedures is expected to **grow** in the coming years.
- RedoTAVR is frequently associated with predicted **CO risk**.
- **In small annuli**, a **SEV** as the index valve is particularly **unfavorable** when the **redo plane is high**. These findings support careful index and redo procedure planning and execution.