

MiRus Siegel™ TAVR

EFS Results

Pradeep Yadav, MD

Presented on behalf of:

Vinod Thourani, Raj Makkar, Aakriti Gupta, Kimberly Holst, Rajiv Gulati, Rahul Sharma, Robert Cubeddu, Dee Dee Wang, Thomas Waggoner, Samir Kapadia, Amar Krishnaswamy, Sahil Khera, Gilbert Tang, Phillippe Genereux, Gennaro Giustino, Martin Leon and Jay Yadav



TCT®

TRANSCATHETER
CARDIOVASCULAR
THERAPEUTICS®

Disclosure of Relevant Financial Relationships

Within the prior 24 months, I have had a financial relationship with a company producing, marketing, selling, re-selling, or distributing healthcare products used by or on patients:

Nature of Financial Relationship

Grant/Research Support

Consultant Fees/Honoraria

Individual Stock(s)/Stock Options

Royalties/Patent Beneficiary

Executive Role/Ownership Interest

Medical Advisory Board

Ineligible Company


Edwards Lifesciences, Abbott, Medtronic, Boston Scientific, Jena Valve, HighLife, Innovalve, Trisol, Laplace, Croi, TriCare, MiRus

Edwards Lifesciences, Abbott, Boston Scientific, Medtronic

Dasi Simulations, Excision Medical

Dasi Simulations, Trisol, TriCare, PiCardia, Excision Medical

Siegel: A Revolutionary THV System



Delivered Crimped On Balloon Through 8 Fr Expandable Sheath

Delivery System




More Flexion

- Small Profile
- 1:1 Torque → Precise commissural alignment
- Crimped on Balloon prior to insertion through sheath



mSheath
8 Fr



Handle

Crimper



Early Feasibility Study


Objective: Safety and efficacy of Siegel TAVR system in patients with symptomatic severe aortic stenosis, all surgical risks

of patients: 30

of sites: Up to 10 sites in the United States

Independent Core Laboratories CT Analysis, Echo Analysis

Investigational Sites



EFS Baseline Clinical Characteristics

Characteristic	n = 30
Age (year)	76 ± 6.21
Male	17/30 (57%)
NYHA Class ≥II	15/30 (100%)
AVA (Vmax) cm ²	0.73
Mean Gradient (mmHg)	40.36
Velocity (m/s)	4.18
DVI	0.23
LVEF (%)	61
STS Score	2.81%
Low	24/30 (80%)
Intermediate/High	6/30 (20%)

Characteristic	n = 30
HTN	20/30 (67%)
DM	7/30 (23%)
CKD	3/30 (10%)
AFib	7/30 (23%)
CAD	20/30 (67%)
LBBB	1/30 (3%)
Permanent Pacemaker	1/30 (3%)
Severe Nickel Allergy	3/30 (10%)

EFS Procedural Results

Procedural	N = 30
General Anesthesia & TEE	2/30 (7%)
Percutaneous Transfemoral Access	30/30 (100%)
Dilation: Pre Post	6/30 (20%) 10/30 (33%)
Aortic 90% / Ventricular 10% Position or Higher	30/30 (100%)
Commissural Alignment	30/30 (100%)
Sheath insertion to valve deployment time	15 minutes

FIH (5) and EFS (30):

Clinical Outcomes at 30 Days N=35

Outcome	Discharge Siegel 26mm & 23mm
All-cause Mortality	0/35 (0%)
Major Stroke	0/35 (0%)
Annular Rupture	0/35 (0%)
Need for 2 nd Valve	0/35 (0%)
Surgery or intervention related to the device, including aortic valve reintervention ^d	0/35 (0%)
Major Vascular Complications ^a	0/35 (0%)
VARC-3 Type 2-4 Bleeding ^c	0/35 (0%)
Moderate or severe aortic regurgitation	0/35 (0%)
Acute Kidney Injury ^b	0/35 (0%)
New pacemaker	2/35 (5.7%)

^a Vascular Complications is defined as major vascular, access-related, or cardiac structural complication according to VARC-3 Guidelines

^b Renal failure is defined as stage III acute kidney injury: Increase in serum creatinine to $\geq 300\%$ (3 x increase compared with baseline) or serum creatinine of ≥ 4 mg/d (≥ 354 μ mol/L) with an acute increase of at least 0.5 mg/dl (44 μ mol/L)



^c VARC-3 Type 2-4 Bleeding: Fatal bleeding OR bleeding in a critical area or organ, such as intracranial, intraspinal, intraocular, or pericardial necessitating pericardiocentesis, or intramuscular with compartment syndrome OR bleeding causing hypovolemic shock or severe hypotension requiring vasopressors or surgery OR overt source of bleeding with drop in hemoglobin of ≥ 5 g/dL or whole blood of packed red blood cells (RBC) transfusion ≥ 4 units

^d Cardiac reintervention includes any intervention that repairs, alters or replaces a previously operated valve OR balloon aortic valvuloplasty OR Surgical aortic valve replacement OR valve in valve

^e CRF adjudicated data

EFS Siegel TAVR Hemodynamics

N=30





*Corelab assessment

EFS Hemodynamics

N=11

23mm Valve





*Corelab assessment

EFS Hemodynamics


N=19

26mm Valve



*Corelab assessment

Siegel TAVR Deployment



No Foreshortening Allows
Consistent Higher
Placement

+

Shorter Balloon in LVOT


+

Less recoil /need for
oversizing




Should Reduce
Conduction System Injury


Case 1



Case 2




Case 3




Case 4




Case 5




Case 6




Case 7




Case 8




Case 9




Case 10




Case 11




Case 12




Case 13




Case 14




Case 15




Case 16




Case 17




Case 18




Case 19




Case 20




Case 21




Case 22




Case 23




Case 24




Case 25




Case 26




Case 27




Case 28



Case 29




Case 30



Siegel TAVR: Ellipticity Ratio


Ellipticity Ratio =
Major axis / Minor axis



Case 1




Case 2




Case 3




Case 4




Case 5




Case 6




Case 7




Case 8



Case 9



Case 10



	Tricuspid	Bicuspid
Siegel	1.02	1.03
Sapien	1.09	1.12
Evolut	1.2	1.3


Case 11



Case 12



Case 13



Case 14




Case 15




Tchetché D., de Biase C. Bicuspid Aortic Valve Anatomy and Relationship with Devices: The BAVARD Multicenter Registry. Circ Cardiovasc Interv. 2019;12:e007107: 1-10.


Case 16




Case 17




Case 18



Case 19



Case 20



Case 21




Case 22




Case 23




Case 24



Case 25




Case 26



Case 27




Case 28



Case 29



Case 30



+ Biocompatibility

Valve	Material
SIEGEL™	Nickel Free
Market Leader	Nickel, Cobalt

+ Pro-healing

Valve	Pro-healing	Surface
SIEGEL™	+++	Nitride Oxide Coating
Market Leader	---	Nickel, Cobalt

∅ Foreshortening

Valve	Foreshortening
SIEGEL™	0%
Market Leader	27%

SIEGEL TAVR



↑ Radial Force

Valve	Radial Force
SIEGEL™	70 N
Market Leader	35 N

+ Tissue

Valve	Tissue
SIEGEL™	Dry Porcine
Market Leader	Dry Bovine

+ Hemodynamics

Benchtop

Valve	EOA	Gradient
26mm SIEGEL™	3.0 cm²	5.4 mmHg
26mm Market Leader	2.4 cm ²	9.2 mmHg

↓ Access Profile

Valve	Sheath Profile
SIEGEL™	8 Fr
Market Leader	14 & 16 Fr

Conclusions – MiRus Siegel Valve

1-year followup of 5 FIH

30-day followup of 30 US Early Feasibility Study patients

The elegant, precise balloon expandable Siegel TAVR showed promising results in the US EFS:

- No Nickel

- 8 Fr access – crimped on balloon

- No foreshortening ---> precise, high placement

- Low pacemaker rate

Hemodynamics are superior and valve is circular – which likely will result in better durability

All risk categories, pivotal randomized control trial is starting soon