

Novel TAVR-AI algorithm for Valve Selection

Femi Philip MD MRCP

Director SH Program

Gabriella S. Philip, John Ko M.D.; James Schipper M.D.; Micheal Chow M.D. and Joseph Huh M.D.



TCT®

TRANSCATHETER
CARDIOVASCULAR
THERAPEUTICS®

Disclosure of Relevant Financial Relationships

I, **Femi Philip, MD** DO NOT have any financial relationships to disclose.

Disclosure of Relevant Financial Relationships

Within the prior 24 months, I have had a financial relationship with a company producing, marketing, selling, re-selling, or distributing healthcare products used by or on patients:

Nature of Financial Relationship

Grant/Research Support

Consultant Fees/Honoraria

Individual Stock(s)/Stock Options

Royalties/Patent Beneficiary

Executive Role/Ownership Interest

Other Financial Benefit

Ineligible Company

Meta

Background

- TAVR valve choice prioritized procedural safety but life-time management is important.
- Clinical need: Balance safety with long-term outcomes – minimize patient prosthesis mismatch (PPM) and TAV-in-TAV feasibility.
- Can AI driven algorithm valve selection balance procedural safety with lifetime management.
- Our approach: AI algorithmic valve selection validated against the MDT team and 30-day echo outcomes.

Study Objectives

Develop and Validate - CT driven TAVR valve selection algorithm using AI

Compare - Algorithm recommendations vs. MDT team

Assess - Concordance, PPM prevention and TAV-in-TAV feasibility

Human-AI Collaborative Framework

- Framework: Human-in-the-loop system
- AI module autonomously processes pre-procedural MDCT
- AI recommendations are reviewed by the MDT
- Modifications are annotated with structured rationale
- Closed-loop retraining pipeline for model refinement
- Outputs are version-controlled and time-stamped

Methods

Study Design:

Prospective cohort of ~200 patients with CT + TTE + procedural data;

Algorithm tested retrospectively on ~700 patients against two independent vendors' stereotyped CT analyses (external validation, confirming alignment);

Prospectively applied to 200 patients and compared to MDT.

Validation Details:

Internal: 5-fold cross-validation on 700 retrospective cases (92% accuracy for PPM prediction, AUC 0.95; 88% for TTF simulation).

External: Vendor benchmarking for CT-derived inputs; Temporal (prospective 200 vs. MDT) for real-world applicability

Methods

- Primary Endpoints

Concordance with the MDT team decisions (Valve type/ access)

- Secondary Endpoints

PPM avoidance (projected EOA used all discordant cases)

TAV-in-TAV feasibility (proprietary algorithm from the training data set)

Methods

Statistical Analysis

Descriptive stats (means \pm SD, frequencies)

Chi-square/Fisher's exact for associations

Cramer's V for effect size

Relative risk reductions with 95% CI

Power: 80% power, $\alpha=0.05$ to detect $\geq 7\%$ absolute PPM reduction

Tools: Python 3.9 (SciPy, StatsModels)

Algorithm Paradigms

Sizing: S3UR by annular area/diameter;
Evolut by perimeter/diameter

Root morphology: classification (Type I, II, III)


PPM thresholds: (EOAi <0.70
moderate, <0.55 severe)

TAV-in-TAV feasibility: VTC, VTSTJ,
VTA using coronary height, sinus, STJ
clearance

Calcium safety: downsizing or overfill
limits if severe

Access integration: iliofemoral
diameter and tortuosity

TAVR Algorithm Learning System: Paradigm Framework



Model Architecture

Model: Feed-forward supervised neural network (MLP)

Inputs: CT + clinical features

Outputs: Valve type, valve size, procedural strategy


Reinforcement Learning Feedback:

Reward = concordance + optimal echo (no PPM, low MG, no AR)


Penalty = mismatch or poor outcomes

Model updates weights to maximize long-term reward


Neural Network + Reinforcement Learning Schematic




TAVI AI: Model 1



ID Type	Label	Value
1	Polygon Annulus Dimensions - Min. Ø	22.8 mm
	Annulus Dimensions - Max. Ø	28.1 mm
	Annulus Dimensions - Avg. Ø	25.4 mm
	Annulus Dimensions - Area derived Ø	25.2 mm
	Annulus Dimensions - Perimeter derived Ø	25.5 mm
	Annulus Dimensions - Area	495.2 mm ²
	Annulus Dimensions - Perimeter	80.2 mm





ID Type	Label	Value
1	Polygon LVOT Diameter - Min. Ø	23.5 mm
	LVOT Diameter - Max. Ø	29.8 mm
	LVOT Diameter - Avg. Ø	26.6 mm
	Annulus Dimensions - Area derived Ø	26.3 mm
	Annulus Dimensions - Perimeter derived Ø	26.6 mm
	Area	541.6 mm ²
	Perimeter	83.6 mm




ID Type	Label	Value
1	Vessel Length Left Coronary Height	8.2 mm
2	Vessel Length Upper LCA Height	16.0 mm



ID Type	Label	Value
1	Vessel Length Right Coronary Height	13.9 mm




ID Type	Label	Value
1	Diameter Sinus of Valsalva Left	31.7 mm
2	Diameter Sinus of Valsalva Non	29.9 mm
3	Diameter Sinus of Valsalva Right	29.2 mm



ID Type	Label	Value
1	Vessel Length Right Coronary Height	13.9 mm
2	Vessel Length Sinus of Valsalva Height	20.5 mm




ID Type	Label	Value
1	Diameter (Min/Max) Sinotubular Junction Diameter - Min. Ø	26.3 mm
	Sinotubular Junction Diameter - Max. Ø	26.8 mm
	Sinotubular Junction Diameter - Avg. Ø	26.5 mm



ID Type	Label	Value
1	Diameter (Min/Max) Ascending Aorta Diameter - Min. Ø	35.9 mm
	Ascending Aorta Diameter - Max. Ø	36.3 mm
	Ascending Aorta Diameter - Avg. Ø	36.1 mm

TAVI AI: DEMO Model:1



Results- Baseline Characteristics

Characteristic	Value
Number of Patients	200
Age (mean \pm SD)	72 \pm 9 years
Female (%)	47
BSA (mean \pm SD)	1.85 \pm 0.25 m ²
STS Risk Score (mean \pm SD)	4.2 \pm 2.1%
Hypertension (%)	82
Diabetes (%)	35
Coronary Artery Disease (%)	48
LVEF (mean \pm SD)	55 \pm 12%
Annular Area (mean \pm SD)	450 \pm 120 mm ²
Annular Perimeter (mean \pm SD)	75 \pm 10 mm
Coronary Heights Left (mean \pm SD)	14 \pm 3 mm
Coronary Heights Right (mean \pm SD)	16 \pm 5 mm
Calcification (Agatston mean \pm SD)	2200 \pm 1500 AU
Mean gradient	42 \pm 8 mm Hg

Results – Baseline Characteristics

Valve types: 60% S3UR, 40% Evolut FX+



Root morphologies: Type I (40%), Type II (25%), Type III (35%)

BSA mean 1.9 m² (range 1.27–2.5)

Coronary heights: LCA mean 13.5 mm; 20% <12 mm

Calcium: Moderate 45%, Severe 30%

Access: >90% transfemoral feasible; 10% borderline



Results: 30-day Outcomes


- Procedural Success: 98%
- All cause mortality: 1.5%
- Disabling stroke: 1.5%
- Major vascular complications: 5%
- Permanent Pacemaker: 8%
- Mean gradient: 9.6 ± 3.8 mmHg

Results

Algorithm vs Valve Team Concordance (N=200)

Primary Endpoint

- Overall concordance 78% Cramer's $V=0.26$, $p<0.001$;
- AI superiority in discordant cases 63.6% ($\chi^2=90.75$, $p=1.63e-21$, $V=0.67$)




Results

PPM

Absolute reduction: 12% (95% CI 6.5-18%);

Relative reduction: 75% (95% CI 50-91.2%, $p=9.23e-10$)


Group	Moderate PPM (%)	Severe PPM (%)	Total PPM (%)
MDT (32/200)	12	4	16
AI (8/200)	4	0	4



Results

TAV-in-TAV

- Overall feasibility: 90%
- Larger-root phenotype: 97.2% feasible
- Smaller-root phenotype: 75.5% feasible
- Most “spacious” group is highest (~100%), and the most “compact” group is lowest (~66.7%).





Limitations

- Single, small study.
- 30-day follow-up included but longer follow-up is needed.
- AI trained on historical data risks overfitting.
- PPM measured and compared to predicted EOA.
- Validation of decision fidelity
- Cost/ Integration barriers not assessed.

Conclusion

- Algorithm shows high concordance with the MDT team
- Reduced the PPM incidence
- Determine TAV-In-TAV feasibility
- Provides standardized, reproducible decision support
- Potential to augment and not replace clinical judgement

Appendix- Model 2.0



Appendix: Algorithm Backend (Pseudocode)

- Inputs: CT annulus, perimeter, LVOT, coronaries, SOV/STJ, BSA, calcium, access
- If short-wide root or low LCA → favor S3UR
- If tubular root + large BSA → favor Evolut
- Check calcium → limit oversizing if S3UR
- Choose size from manufacturer sizing tables
- Predict EOAI → classify PPM risk
- Output: Valve system, size, procedural notes
- Feedback: Compare to Valve Team choice + 30-day echo
- Reward = concordance + optimal hemodynamics
- Reinforcement loop updates weights to favor high reward patterns