

Technical Considerations for TAVR in Bicuspid AS

Giuseppe Tarantini, MD, PhD

Director of Interventional Cardiology Unit, University of Padova



TCT[®]

TRANSCATHETER
CARDIOVASCULAR
THERAPEUTICS[®]

Disclosure of Relevant Financial Relationships

Within the prior 24 months, I have had a financial relationship with a company producing, marketing, selling, re-selling, or distributing healthcare products used by or on patients:

Nature of Financial Relationship

Consultant Fees/Honoraria

Companies

Abbott Laboratories, Boston
Scientific, Edwards Lifesciences,
Medtronic, GADA, Microport, SMT

OPERATIVE OVERVIEW

☐ PHENOTYPING

☐ SIZING

☐ POSITIONING

☐ OPTIMIZATION

OPERATIVE OVERVIEW

☒ PHENOTYPING


☐ SIZING

☐ POSITIONING


☐ OPTIMIZATION

BAV MORPHOLOGY CHARACTERIZATION


SIEVERS (SURGICAL-derived)¹




Bicommissural
NON raphe-type



Bicommissural
raphe-type







Tricommissural
with calcific fusion



JILAIHAWY (Operative CT-derived)²

CT-SCAN



OPERATIVE OVERVIEW

☐ PHENOTYPING

☒ SIZING


☐ POSITIONING

☐ OPTIMIZATION

CT SCAN based #1


VIRTUAL BASAL RING (Annular)

***Challenging in Type 0 BAV**




VIRTUAL RAPHE RING (Supra-annular)

***Multi-parametric evaluation**





VBR vs. VRR sizing




VIRTUAL RAPHE RING		Value
Perimeter (mm)		78.7
Intercommissural length (mm)		28.0
Ante-Raphe length		20.5
Perimeter-derived diam. (mm)		25.1
Perimeter (mm)		78.7
Intercommissural length (mm)		28.0
VIRTUAL BASAL RING		Value
Perimeter (mm)		89.8
Min Diam (mm)		27.2
Max Diam (mm)		28.4
Mean Diam (mm)		27.8
Perimeter-derived diam. (mm)		28.6
Area (mm ²)		628

CT SCAN based #2


Type 0 Tricommissural



CO-DOMINANT




ANNULAR-DOMINANT




VBR sizing


Type 1 and 2




CO-DOMINANT



ANNULAR-DOMINANT



RAPHE-DOMINANT




VBR sizing

***VRR
sizing***

BAV based #3

Balloon valvuloplasty to test the supra-annular structure

- ❑ **Supra-annular vs Annular WAISTING**
- ❑ **Symmetrical vs. Asymmetrical EXPANSION**







THV TYPE CHOICE



Evolut R/Pro/Pro+/Fx


CE mark for BAV from June 2020

- 
- Repositionable/Retrievable
 - Lower risk of annular rupture^{1,2}
 - Supra-annular valve circular
 - Better hemodynamics^{1,2}

- 
- More PVL^{1,2}
 - Higher risk of Pacemaker^{1,2}
 - Impaired Coronary Access

Sapien 3/Ultra/Resilia

Caution for BAV removed from IFU

- 
- Higher radial force
 - Maintains circular shape
 - Lower PVL^{1,2}
 - Coronary «friendly»

- 
- Higher risk of annular rupture^{1,2}
 - Higher risk of aortic dissection

OPERATIVE OVERVIEW










☐ PHENOTYPING

☐ SIZING


☒ POSITIONING

☐ OPTIMIZATION

Annular- Vs. Supra-annular Positioning




Supra-annular Positioning shades



BEV

- ❑ **THV embolization risk**
- ❑ **Prosthesis-patient mismatch (down-sizing)**



SEV

- ❑ **Coronary access impairment²**
- ❑ **Redo-TAVR feasibility¹**

OPERATIVE OVERVIEW

☐ PHENOTYPING

☐ SIZING

☐ POSITIONING

☒ OPTIMIZATION

THV expansion assessment

FLUOROSCOPIC

EVOLUT PRO THV

ACURATE Neo THV







PORTICO THV

LAO/CRAN

Co-planar view

RAO/CAU


Orthogonal view




INVASIVE

**EVOLUT 29-mm
RAPHE-TYPE BAV
P-P Grad 40 mmHg**

**Post-dilatation
23-mm CB
P-P Grad 10 mmHg**



Post-dilatation with SEVs



Take home message

BAV is NOT a MONOLITIC ENTITY

- ❑ TAVR for **NON-RAPHE TYPE** and **FUNCTIONAL BAV** resemble classic TAV
- ❑ **TAILORED** THV SIZING and POSITIONING is the **RULE** in **RAPHE-TYPE**
- ❑ **POST IMPLANT ASSESSMENT AND OPTIMIZATION** IS RECOMMENDED



@interventionalcardiologyteampadova



@giuseppetarantinimdphd



@G_Tarantini01