

# Dedicated Technologies for Treatment of Aortic Regurgitation (AR)

Nick Amoroso, MD FSCAI FACC

Medical University South Carolina

---



**TCT**®

TRANSCATHETER  
CARDIOVASCULAR  
THERAPEUTICS®

# Disclosure of Relevant Financial Relationships

Within the prior 24 months, I have had a financial relationship with a company producing, marketing, selling, re-selling, or distributing healthcare products used by or on patients:

## Nature of Financial Relationship

Grant/Research Support

Consultant Fees/Honoraria

Individual Stock(s)/Stock Options

Royalties/Patent Beneficiary

Executive Role/Ownership Interest

Other Financial Benefit

## Ineligible Company

Edwards Lifesciences, JenaValve, Vdyne, Tricare, P&F

Abbott, Boston Scientific, Edwards Lifesciences, JenaValve, Vdyne, Egnite, EastEnd Medical

Nininger Medical

n/a

n/a


n/a

# Aortic Regurgitation – The Problem


- Annulus size often large



- Younger patients than aortic stenosis; Lifetime management (valve & aneurysms)



# Preclinical Device - Cusper






# TAVR Devices

	J-Valve	Trilogy
Manufacturer	JC Medical	JenaValve Technology
Approval	NMPA (2017)	CE (2021)
Mechanism of expansion	Self-expanding	Self-expanding
Leaflet material	Bovine pericardial	Porcine pericardial
Leaflet position	Intra-annular	Supra-annular
Frame material	Nitinol	Nitinol
Frame height	17–25 mm*	32–36 mm*
Frame cell size	Sinus cut-outs	27–31 Fr
Native leaflet interaction	U-shaped anchor rings	Locator
Commissural alignment	Self-aligning design	Self-aligning design
Sealing	Fabric, anchor rings	Flared sealing ring, locators
Access	Transfemoral	Transfemoral
Sizes	22, 25, 28, 31 and 34 mm	23, 25 and 27 mm
Target annular diameter range	18.0–33.1 mm	21.0–28.6 mm
Target annular perimeter range	57–104 mm	66–90 mm
Delivery system flexibility/steerability	+/++	++/++
Repositionable system	+	+
Retrievable system	–	–
Delivery sheath size compatibility	18–22 Fr (ID) or 16 Fr Edwards ESheath <sup>†</sup>	Dedicated 20 Fr (ID) <sup>‡</sup>

## JenaValve Trilogy TAVR Device



Locator grasping of leaflets



Valve Deployment

## J-Valve TAVR Device



U-shaped anchor rings positioned on the leaflets



Valve Deployment

to  
3))

# Meta-Analysis of Dedicated Versus Off-Label TAVR for Native Aortic Valve Regurgitation

aya et al.  
2017<sup>20</sup>  
41  
72  
42

n  
Age, year  
Female, %

## Challenges of TAVR in native AR

## TAVR and THV design specific solutions

Larger valve sizes

**Table 2** Subgroup analysis of clinical outcomes in on-label devices, off-label SE and off-label BE: a meta-analysis

	On-label devices			Off-label SE			Off-label BE				
Outcomes	Event rate (95% CI)	<i>I</i> <sup>2</sup> (%)	P value*	Event rate (95% CI)	<i>I</i> <sup>2</sup> (%)	P value*	Event rate (95% CI)	<i>I</i> <sup>2</sup> (%)	P value*	$\chi^2$	P value†
In-hospital Outcomes											
All cause mortality	0.02 (0.01 to 0.03)	4.7	0.398	0.04 (0.02 to 0.08)	11	0.346	0.04 (0.02 to 0.07)	0	0.797	5.52	0.063
Technical success	0.97 (0.94 to 0.98)	11.8	0.329	0.85 (0.78 to 0.90)	26.7	0.234	0.92 (0.88 to 0.95)	21.4	0.280	17.46	0.000
Device success	0.95 (0.92 to 0.97)	0	0.598	0.83 (0.76 to 0.87)	66.6	0.001	0.89 (0.76 to 0.96)	88.9	0.003	19.75	0.000
PPI	0.10 (0.06 to 0.15)	74.3	0.000	0.19 (0.14 to 0.24)	34.3	0.133	0.18 (0.10 to 0.29)	80.7	0.006	6.18	0.046
Moderate or severe AR	0.02 (0.01 to 0.03)	0	0.745	0.04 (0.02 to 0.07)	31.6	0.156	0.08 (0.06 to 0.12)	1.1	0.364	27.05	0.000‡§
SVI	0.02 (0.01 to 0.03)	18.4	0.269	0.15 (0.11 to 0.21)	46.2	0.053	0.05 (0.03 to 0.08)	0	0.382	50.92	0.000‡§¶
Valve migration	0.02 (0.01 to 0.04)	10.8	0.346	0.10 (0.05 to 0.22)	84.7	0.000	0.07 (0.04 to 0.11)	42.1	0.178	10.89	0.004‡§
Major bleeding	0.04 (0.02 to 0.06)	0.2	0.439	0.05 (0.03 to 0.09)	74.2	0.000	0.03 (0.01 to 0.08)	76.9	0.005	0.86	0.651
AKI 2 or 3	0.02 (0.01 to 0.04)	0	0.647	0.04 (0.01 to 0.17)	50.1	0.138	0.03 (0.01 to 0.18)	81.9	0.019	57.9	0.006
30-days Outcomes											
All cause mortality	0.03 (0.02 to 0.05)	41.7	0.080	0.06 (0.03 to 0.11)	20.3	0.285	0.06 (0.03 to 0.11)	0	0.537	4.14	0.162
Stroke	0.02 (0.01 to 0.05)	0	0.938	0.05 (0.03 to 0.07)	0	0.938	–	–	–	2.07	0.151
Moderate or severe AR	0.01 (0.00 to 0.03)	0	0.870	0.09 (0.03 to 0.23)	75.1	0.007	–	–	–	8.02	0.005

Regurgitation. Interventional Cardiology 2025;20:e25.

30-day mortality

Log EuroSCORE = 24.2% ± 0.3%


Peng Y, Lin Y, Su Z, et al. Comparative outcomes of on-label and off-label transcatheter aortic valve replacement for aortic regurgitation: a systematic review and meta-analysis. Open Heart 2025;12:e003482.

■ Dedicated Devices ■ Off-Label Devices

CRF®  
TCT

# J-Valve Early Feasibility Study – 15 patients

- No procedural mortality, coronary obstruction, embolization, migration, or ViV
- One convert to surgery for failed device release with tortuous aorta
- One 30-day noncardiac death
- LV remodeling at 30-days



**TABLE 1** Echocardiographic Characteristics

	Full Cohort Baseline (n = 15)	Paired Cohort <sup>a</sup>		P Value <sup>b</sup>
		Baseline (n = 14)	30 d (n = 14)	
LVEF, %	53.84 ± 7.97	54.10 ± 8.20	49.06 ± 9.31	0.064
AV mean gradient, mm Hg	5.38 ± 2.20	5.13 ± 2.05	5.57 ± 2.04	0.625
EOA, cm <sup>2</sup>	3.04 ± 0.68	3.11 ± 0.65	2.90 ± 0.68	0.444
Residual AR severity ≥mild	NA	NA	0	NA
Paravalvular regurgitation	NA	NA	0	NA
LVEDD, cm	6.00 (5.10-6.70)	5.90 (5.10-6.20)	5.20 (4.80-5.50)	0.014
LVESD, cm	4.20 (4.00-5.60)	4.20 (4.00-5.20)	3.95 (3.00-4.50)	0.088
LVEDV, mL	167.70 (131.80-232.10)	166.95 (131.80-229.00)	133.10 (109.50-201.10)	0.017
LVESV, mL	87.70 (53.70-115.90)	80.00 (53.70-105.50)	65.60 (51.30-116.50)	0.241 <sup>c</sup>
LV mass, g	222.00 (157.00-287.00)	219.50 (157.00-272.00)	189.00 (165.00-236.00)	0.056

# JenaValve ALIGN AR

## Transcatheter

### Death

#### Any stroke

Disabling stroke

Non-disabling stroke

#### Major or life-threatening bleeding

#### Major vascular complication

#### Acute kidney injury stage 2 or 3 or dialysis (7 days)

#### Surgery or intervention related to the device

Aortic Endograft and Commercial THV for aortic dissection

Surgical aortic valve replacement for Trilogi transcatheter

Commercial transcatheter heart valve for Trilogi transcatheter

Trilogi transcatheter heart valve for Trilogi transcatheter

#### New pacemaker implantation

#### Moderate or greater aortic regurgitation

#### Total

Data are n (%) or n/N (%). \*p < 0.05; †p < 0.01; ‡p < 0.001; §p < 0.0001

Table 2: Primary safety endpoints at 30 days

Achieved the  
(7.8% vs. 25

346 patient

### ALIGN-AR 2024<sup>40</sup>

n	180
Age, years	75
Female, %	47
STS-PROM, %	4

#### Trilogi

#### THV

Procedural/ device success	95
----------------------------------	----

Residual AR >moderate	6
--------------------------	---

THV migration/ embolisation	2
--------------------------------	---

Second THV	2
------------	---

Surgical conversion	1
------------------------	---

New PPM	24
---------	----

Stroke/TIA	2
------------	---

Major vascular complication	4
--------------------------------	---

Major bleeding	4
----------------	---

30-day mortality	2
---------------------	---

anatomic criteria

regurgitation severity grade 3 or

firm

diameter <66 mm or >90 mm

or

aorta

regurgitation greater than

moderate aortic stenosis (mean

gradient >50 mm

Hg)

or

coronary length <55 mm

converted to surgical aortic valve

replacement

received commercial aortic valve

replacement

180 enrolled


177 Trilogi

valve st

177 included in 1-year valve implantation


study population

Class III Class IV



Timepoint

### C 6-min walk test






# JenaValve – ALIGN AR Trial + CAP



## Screening and Patient Disposition



\* Available as of 3/22/2025

# Take-Home Message

- Continued Unmet Needs for Treatment of Aortic Regurgitation
- Off-label TAVR Suffers Unacceptable Outcomes at Scale
- Dedicated Transcatheter Therapies (e.g. JenaValve, J-Valve)
  - Safer
  - More technical success
  - Less PVL or short-term reintervention
- **Needs:** Commercial Availability, Pacemaker Avoidance, More Anatomies, Trials of all risk profile, Trials vs surgery