

Transcatheter Aortic Valve Replacement with Ken-Valve in Aortic Regurgitation (Jenscare Scientific): One Year Outcomes And Challenging Cases

Anson Cheung, MD

St. Paul's Hospital, Vancouver, Canada

On-behalf of Clinical Trial Team



TCT[®]

TRANSCATHETER
CARDIOVASCULAR
THERAPEUTICS[®]

Disclosure of Relevant Financial Relationships

Within the prior 24 months, I have had a relevant financial relationship(s) with an ineligible company(ies) listed below.

<u>NATURE OF FINANCIAL RELATIONSHIP</u>	<u>INELIGIBLE COMPANY</u>
Grant/Research Support	Edwards Lifesciences, Shockwave, Abbott Vascular, Corysm
Consultant Fees/Honoraria	Medtronic, Edwards Lifesciences, Abbott Vascular, Boston Scientific, Jenseare
Individual Stock(s)/Stock Options	Total Flow Medical Shockwave Kardium Huihe Vesalius
Other	Eligibility Committee Board Member (TRINITY Trial)

Innovative Design of Ken-Valve TAVR Device for AR



Precise deployment & stable anchoring

Integrated clamps & radial force

- Stronger leaflet clamping for stable anchoring
- Enhanced radial force to reduce valve migration
- Minimal oversizing required



Clamps wrapped by bovine pericardium to minimize tissue injury



Anchor marker ensures precise positioning




Soft PTFE skirt to prevent PVL and to reduce conduction defect

Supra-annulus leaflet design for optimal EORA



Short stent frame for future coronary access

Ken-Valve TAVR Delivery System




Ken-Valve TAVR Procedure



MKTMM2025051501

Ken Valve Pivotal Trial



Outcome Measure : MVARC criteria

- Study Design

- **Prospective**
- **Single Arm**
- **Multi-center Study in China** (15 enrolling sites)

- Objective

To evaluate the safety and efficacy of the Ken-Valve TAVR System in treating symptomatic **AR** (or combined with stenosis) patients with high surgical risk.

- Endpoint

All cause mortality, heart function, quality of life, device success rate at 12 months, and MACE event at 30 days.

- **Follow up:** Discharge, 30d, 6m, 12m

■ Coordinating Investigator

The First Affiliated Hospital of Naval Medical University | Prof. XU Zhiyun
West China Hospital, Sichuan University | Prof. GUO Yingqiang



Baseline Characteristics

Clinical Characteristics (n=142)

Age — yr	70.3 ± 5.5	COPD	40.1% (57/142)
Male Gender	63.68% (90/142)	Peripheral vascular disease	48.6% (69/142)
BMI — kg/m2	23.1 ± 3.3	Prior surgical history	51.4% (73/142)
Moderate to severe frailty	62% (88/142)	NYHA class	
Diabetes	11.3% (16/142)	I	0
Hypertension history	71.8% (102/142)	II	0.7% (1/142)
Coronary artery disease	43% (61/142)	III	47.9% (68/142)
Previous stroke or TIA	21.8% (31/142)	IV	51.4% (73/142)
Atrial fibrillation	11.3% (16/142)	STS score	5.9 ± 3

Distribution of Ken-Valve Sizes Used

Ken-Valve Sizes Used In The Trial



✓ **46% of enrollments utilized valve sizes > 29mm**

Baseline Characteristics and Procedural Outcomes

Echocardiographic Data (N=142)		Procedural Outcomes (N=142)	
Severity of aortic regurgitation		Procedural time (mins)	121 ± 63
severe, grade 3+	100%	Device time (mins)	9 ± 9
Severity of aortic stenosis		Procedural Success	96.5%**
mild, grade 1+	9.2%	Device Success	97.2%*
moderate, grade 2+	2.1%	Conversion to surgery	1.4% (2/142)
severe, grade 3+	12%	Acute renal dysfunction	0.7% (1/142)
Effective regurgitant orifice area (cm ²)	0.5 ± 0.3	Needed for CPB	0.7% (1/142)
Vena contra width (cm)	0.8 ± 0.2	In-hospital or 30-day mortality	2.1%(3/142)
Stroke volume (ml)	77 ± 35		
Left ventricular end-systolic dimension (cm)	45 ± 11		
Left ventricular end-diastolic dimension (cm)	62 ± 10		
Average transvalvular gradient (mmHg)	33 ± 18		
Peak transvalvular velocity (m/s)	3.6 ± 1.1		
Left ventricular ejection fraction mean (%)	54 ± 13		


* Device Success: 4 patients needed a second valve or open heart surgery due to misposition of the Ken-Valve.

Major Adverse Events at 1 Year

All-cause mortality	5.6% (8/142)
Permanent pacemaker implantation rate	14.1% (20/142)
III degree AV block	8.5% (12/142)
Other types of arrhythmias	5.6% (8/142)
Major bleeding	4.2% (6/142)
Stroke	2.1% (3/142)
Reintervention	0.7% (1/142)
Coronary artery obstruction	0
Thrombosis	0
Cardiac tamponade	0


Significant AR Reduction at 1 Year

94.7% of the patients had mild or less AR or PVL at 12 months ($P < 0.001$)




Cardiac Function and Quality of Life at 1 Year

NYHA Classification




96% of the patient in NYHA Class I/II at 12 months ($P<0.001$)

EQ-5D



Significant improvement in QoL at 12 months ($P<0.001$)

EOA and Mean Transvalvular Gradient at 1 year



Case I

- 67 y/o Male in NYHA Class III, diabetes, post PCI
- Enlarged left ventricle with LVEF - 28%, PASP=50mmHg
- Severe aortic valve regurgitation (4+), no mitral or tricuspid valve disease


Case I



Key challenges of the case:

1. Severely dilated annulus and STJ (mean annular diameter - 30.5mm and 47.4mm)
2. Poor LVEF - 28%.



Case 1



Procedural steps:

1. Advancement of the delivery system over a guidewire, flexion of the delivery system to achieve co-axiality
2. Deployment of the three clamps into aortic sinus
3. Valve deployment
4. Stable valve position with trivial PVL

1 year follow up




1. NYHA I/II
2. Well functioning Ken-Valve with excellent hemodynamics and no PVL
3. LVEF improvement from 28% to 51%
4. LV remodeling from 229/166ml to 111/54ml

Case II

- 80 y/o male in NYHA Class III
- Hx of hypertension, pneumonia, post and prostate and esophagus cancer surgery and radiotherapy
- Enlarged left ventricle and left atrium, LVEF=44%
- Severe aortic valve regurgitation (4+), mild to moderate mitral regurgitation


Case II



Key challenges of the case:


1. Dilated aortic annulus and STJ (mean diameter - 30.3mm and 41.9mm)
2. Horizontal heart, co-axiality could be difficult to achieve

Case II



Follow-up at 1 month:

1. Well functioning well valve with no PVL
2. Clinical improvement in NYHA class and QoL



Conclusions

- **Ken-Valve TAVR system is an innovative device, designed specifically to treat aortic regurgitation (AR) or predominantly AR mixed with AS**
- **Device advantages:**
 - ✓ Integrated clamping + radial force anchoring mechanism, provide reliable anchoring
 - ✓ PVL mitigating valve stent design
 - ✓ Multiple and large valve size options to accommodate various anatomy
 - ✓ Deflectable delivery system allows better co-axiality and alignment

Conclusions

- **Ken-Valve Pivotal Trial demonstrated good procedural safety and effectiveness**
- **Excellent device performance with large EOA, low transvalvular gradient with minimal PVL**
- **Significant clinical improvement in NYHA class and QoL with low rate of adverse events**
- **Encouraging treatment option for high-risk patient with severe AR and mixed AR**



THANK YOU



TCT®

TRANSCATHETER
CARDIOVASCULAR
THERAPEUTICS®